

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-004885, in compliance with the agency's Maintenance Policy and Procedures.

Signed

Joshua Smith

Name of Position: Director/Records Custodian

JORDAN FAITH THOMPSON NOTARY PUBLIC STATE OF OKLAHOMA COMMISSION NO. 25010889 EXPIRES 09-16-2029 State of Oklahoma

County of Oklahoma

Signed or attested before me on 10/00/2025

by Joshua Smith.

(Signature of notarial officer)

Service Overview

3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 4/14/2025

SERIAL: 80-004885

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- BOT Dry Gas Certificate of Analysis
- Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis on Certified Solutions (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seq) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- Maintenance was performed on the Intoxilyzer 8000, Serial No. 80-004885 on 4/14/2025 in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 et seq) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
 - Manufacturer: ILMO Unit: 208 ppm / 105 L
 - Cylinder Lot No. Canister No.: 00925080A4-018
 - BOT analysis average test result: 0.082
 Cylinder expiration date: 2/5/2027

3. In accordance with 40 O.A.C. 25-1-2.1, the approved reference standards and equipment were used:

Simulator Model			
Guth 2100	Guth 2100	Guth 2100	Guth 2100
Serial #		had fight and	
DR3591	DR3694	DR2455	DR3641
Concentration			
.02	.04	.10	.20
LOT#			
24150	24330	24390	24080
Manufactured Date	A Atlanta	4 - 44	
4/16/2024	9/10/2024	10/29/2024	2/13/2024
Expiration Date			
4/16/2026	9/10/2026	10/29/2026	2/13/2026
Solution Commission Date			· · · · · · · · · · · · · · · · · · ·
4/7/2025	4/14/2025	4/14/2025	4/14/2025

Completed by BOT Personnel:

Christopher Pape
Signature/Name

Breath-Alcohol Testing Program Administrator

Title

Certified by:

Signature

Joshua Smith Name State Director Title

ENCH CHECK DATE 04/14/2025	TIME PROCEDURE BEGAN		INSTRUMENT LOCATION 80-004885	5	
STRUMENT TYPE	INSTRUMENT		201 011 212 21	INSTRUMENT SE	
Intoxilyzer	I-80		TOUR TO	80-004	1885
		IESI KI	ESULTS		
Diagnosti PASS					
Mouth Alcoh	ol Test				
PASS					
RFI Tes					
PASS					
Abort Te		Di	uplicate Breath	Sampling Re	sults
PASS		Value	Test 1	Test 2	Control Test
Improper Sam	\$1000 H 360 W 1		1 2 2 2 2		0 017
PASS	743.2	.02	0.017	0.018	0.017
Interferent Det		.04	0.038	0.037	0.037
PASS		.10	0.096	0.094	0.096
No Sample Giv		.20	0.193	0.192	0.193
PASS		.20			0.155
Insufficient San PASS	1 (1) (1) (1) (1) (1) (1)		PA	SS	
Air Blank Contam PASS					
Range Exceed					
		CERTIFYING	TECHNICIAN		
I swear under penalt	ty of perjury that in	regards to th	ne above listed	instrument, l	have complied
the Intoxliyzer Bend					
GNATURE				DATE	
MATE .	Bn	-		4	14/20)
PAPE, CHRISTO	PHER			PERMIT#	0023020



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT	1000	
Intoxilyzer - Alcohol	Analyzer	
Model 8000	SN	80-002591
04/10/2025		10:17

Test	g/210L	Time
Air Blank	0.000	10:17
Cal Check	0.082	10:18
Air Blank	0.000	10:18
Cal Check	0.082	10:19
Air Blank	0.000	10:19
Cal Check	0.082	10:19
Air Blank	0.000	10:20
Cal Check	0.082	10:20
Air Blank	0.000	10:21
Cal Check	0.082	10:21
Air Blank	0.000	10:22
Cal Check	0.082	10:22
Air Blank	0.000	10:22
Cal Check	0.082	10:23
Air Blank	0.000	10:23
Cal Check	0.082	10:24
Air Blank	0.000	10:24
Cal Check	0.082	10:24
Air Blank	0.000	10:25
Cal Check	0.082	10:25
Air Blank	0.000	10:26
Cal. Check Stats		
Average	0.0820	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

THE RESERVE OF THE PARTY OF THE	0.08 g/210L = 208 PPM ÷ 2605*	Log to obe
Manufacturer: ILMO	Manufacturer Cert. No.: 17033	Cylinder Size: 105L
Lot No.: 00925080A4	Cylinder No.: 018	Expiration Date: 02/05/2027
calibrating unit in accordance with 47 O.S.	. 759, 40 OAC 25-1-3, and appears on the current confe	reference standard meets the requirements for use as a
Anal. Toxicol. (1996) 20(6): pp. 484-491.	of Transportation in the Federal Register [77 FR 64588	et seq]. * See Dubowski, K.M. and Essary N., Journal of
Anal. Toxicol. (1996) 20(6): pp. 484-491. Signature	of Transportation in the Federal Register [77 FR 64588	et seq]. * See Dubowski, K.M. and Essary N., <i>Journal of</i>



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID:

17033

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

00925080A4

Expiration:

2/5/2027

RCVD 2/13/25

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Accuracy

Analytical

Component:

Concentration:

(U, k=2): abs

Method:

Ethanol Nitrogen 208 ppm Balance +/-0.002 BAC(G/210L) NDIR

[5.2 ppm]

Distributed by:

CMI Inc.

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

*Reference Standard:

NIST Traceable Certified Reference Material - 262.3 µmol/mol Ethanol in Nitrogen - Serial No. ND7017 Lot No. 080722E2

Store in dry area, away from sources of heat. ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

01-22-2025

Issuance Date



The calibration results within this certificate were obtained at the facility listed above using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service. This certificate applies only to the items described and shall not be reproduced other than in full, without written approval from the issuing facility.

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

Section 1 ARREST TIME CITATION# ARREST DATE COUNTY # CITY# 04/14/2025 1500 55 COUNTY ARREST LOCATION BOT OKLAHOMA DATE OF BIRTH HEIGHT SUBJECT NAME 08/09/1961 510 GAYLON LEVERNE M 165 NORTH, ZIP CODE ADDRESS 73133 OKLAHOMA CITY 3600 MARTIN LUTHER KING A OK DRIVER LICENSE # XPIRATION DATE CLASS 06/30/2019 T083454871 OK A M VEHICLE MAKE PLACARD MOCK MOCK MOCK OK N N

On the above date, time, and location, the above named person was arrested, and i had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, tumpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law.

(Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE, INSTALLED NEW CYLINDER

THE PERSON WAS INFORMED OF THE IMPLIED CONSENT TEST REQUEST BREATH-ALCOHOL ANALYSIS RECORD, REPORT OF THE PERSON NAMED IN SECTION 1 ABOVE, and LOG OF TEST AND MAINTENANCE RECORD. Section 2 SERIAL NUMBER SPECIALIST I-8000 BOARD OF TESTS 80-004885 STANDARD ARGET VALUE PERMIT NUMBER Dry MANUFACTURER 0.080 23020 DATE CYLINDER INSTALLED DATE CYLINDER EXPIRES 04/14/2025 00925080A4 02/05/2027 ILMO OBSERVATION START TEST DATE A 15-minute (or greater) deprivation/observation 04/14/2025 1520 period was observed prior to and continuing RESULT TABLE to the time of test administration. Time Test g/210L I administered said Breath-Alcohol Test in 15:48 0.00 Air Blank Subject Test 0.00 15:49 compliance with the applicable operating Breath Volume 2.066 LITERS procedure of the Board of Tests. Air Blank 0.00 15:49 15:51 Wait Air Blank 0.00 15:51 Subject Test 0.00 15:52 SIGNATURE OF OPERA 1.972 LITERS Breath Volume 0.00 15:52 Test Result Cal Check 0.07 15:53 PAPE, CHRISTOPHER Air Blank 0.00 15:53 BrAC g/210L 0000023020 0000000484 0.00 AGENCY BOARD OF TESTS

What happens next?

Section 3

- Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges.
- 2. If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle.

 3. You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
- 4. The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma.

To confirm or change your address on file, contact Service Oklahoma.

5. The Order and Notice of Revocation will inform you of your requirements to regain driving privileges in Oklahoma.

In accorda	ince with Title 1	2 O.S. Section 426, "I state under	penalty of perj	ury under the laws of Oklahoma that the for	regoing is true and correct."	Section 4
Date:		Place (location when signed):		Signature of arresting office	r	الششاشا
PAPE	, CHRIS	TOPHER		BOARD OF TESTS		
BADGE 0000	000484	AGENCY ADDRESS INTER-AGENCY		ZIP CODE 0 0 0 0 0	4054252460	
OTHER W	VITNESSES:					Section 5
1 NAME	В		TITLE	ADDRESS	PHONE	
2 NAME	E		TITLE	ADDRESS	PHONE	

8/17

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

							Section
ARREST DATE 04/14/2025	1500		55 ci	IY#	CITATION		
ARREST LOCATION BOT			CITY			OKLAH	OMA
SUBJECT NAME NORTH, GAYLON LEVE	RNE			09/19	61 M	неідн 51	12 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2
ADDRESS 3600 MARTIN LUTHER	KING A	OKLAHO	MA CITY		OK	73	DE 133
T083454871			730/2019	OK OK	A	RESTRICTIONS	M
VEHICLE MAKE MOCK	MOCE	ζ		MOCK	7	OK N	N N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, tumpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law.

(Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE, INSTALLED NEW CYLINDER

		THE PERSON V	VAS INFORMED OF	THE IMPLIED CO	DISENT TES	ST REQUEST	
BREATH-ALCOHOL ANALYSIS	S RECORD, REPORT OF THE F	PERSON NAMED	IN SECTION 1 ABOV	E, and LOG OF T	EST AND M	AINTENANCE RECORD.	Section 2
I-8000			NUMBER - 004885	SPECIALIST BOARD	OF T	TESTS	
STANDARD Dry			0.080		PERMIT N		
MANUFÁCTURER ILIMO			00925	080A4	DATE CYL	LINDER INSTALLED	DATE CYLINDER EXPIRES 02/05/2027
	T DATE 4/14/2025					ation/observation	
RESULT TABLE Test	g/210L	Time		the time of to		and continuing nistration.	OF THE STATE
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume Air Blank Cal Check Air Blank	0.00 0.00 2.066 LITERS 0.00 0.00 0.00 1.972 LITERS 0.00 0.07	15:48 15:49 15:49 15:51 15:51 15:52 15:52 15:53 15:53	complian proc signature signature NAME PAPE, CHI BADGE #	edure of the	e applicate Boar		Test Result BrAC g/210L
			AGENCY BOARD O	L. P. Brander		0023020	0.00

What happens next?

Section 3

- Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges.
- 2. If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle.
- 3. You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
- The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma.
 To confirm or change your address on file, contact Service Oklahoma.
- 5. The Order and Notice of Revocation will inform you of your requirements to regain driving privileges in Oklahoma.

In accordance with Title 12 O.S. Section 426, "I state under	r penalty of perjury under the laws of Oklahoma that the foregoing is true and correct."	Section 4
Date:Place (location when signed):_	Signature of arresting officer:	
PAPE, CHRISTOPHER	BOARD OF TESTS	
BADGE AGENCY ADDRESS INTER-AGENCY	ZIP CODE	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma Cily, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT		
Intoxilyzer - Alcohol	Analyzer	
Model 8000	SN	80-002591
01/15/2025		09:13

Test	g/210L		Time
Air Blank	0.000		09:14
Cal Check	0.021		09:14
Air Blank	0.000		09:15
Cal Check	0.021		09:16
Air Blank	0.000		09:16
Cal Check	0.021		09:17
Air Blank	0.000		09:17
Cal Check	0.021		09:18
Air Blank	0.000		09:19
Cal Check	0.020		09:19
Air Blank	0.000		09:20
Cal Check	0.021		09:21
Air Blank	0.000		09:21
Cal Check	0.021		09:22
Air Blank	0.000		09:22
Cal Check	0.021		09:23
Air Blank	0.000		09:24
Cal Check	0.021		09:24
Air Blank	0.000		09:25
Cal Check	0.020		09:26
Air Blank	0.000		09:26
Cal Check Stats			
Average	0.0208		
Std Dev	0.0004	4	
Rel Std Dev(%)	2.0270		

Manufacturer: Guth	Manufactured Date: 4/16/2024	Concentration: 0.020
Lot No.: 24150	Expiration Date: 4/16/2026	Equipment Model: Guth 2100
prescribed by the State Director of Tests ports 40 OAC 25-1-2.1. Approved breath-alcol	Infrared Spectroscopy utilizing the Intoxilyzer 8000 and ursuant to 40 OAC 30-1-3. The above NIST traceable rand measurement equipment and reference solutions. E the for Breath Alcohol Testers (72 FR 34747), published	eference solution meets the requirements for use pursuan quipment used appears on the current or supplemented



RCVD 12/3/24 CD

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 14, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2419% (w/vol) ethyl alcohol. The expiration date for this lot number is February 13, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03132302 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

10:32

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-002591 04/07/2025

Test	g/210L	Time
Air Blank	-0.000	10:32
Cal Check	0.040	10:33
Air Blank	0.000	10:34
Cal Check	0.040	10:34
Air Blank	0.000	10:35
Cal Check	0.040	10:36
Air Blank	0.000	10:36
Cal Check	0.040	10:37
Air Blank	0.000	10:37
Cal Check	0.040	10:38
Air Blank	0.000	10:39
Cal Check	0.040	10:39
Air Blank	0.000	10:40
Cal Check	0.040	10:41
Air Blank	0.000	10:41
Cal Check	0.040	10:42
Air Blank	0.000	10:42
Cal Check	0.039	10:43
Air Blank	0.000	10:44
Cal Check	0.039	10:44
Air Blank	0.000	10:45
Cal Check Stats		
Average	0.0398	
Std Dev	0.0004	
Rel Std Dev(%)	1.0593	

NIST Traceable Breath-Alcoho	l Equivalent (g/210L +/- 3%) Ratio of ethanol and w	vater free of test interfering substances.	
Manufacturer: Guth	Manufactured Date: 9/10/2024	Concentration: 0.040	
Lot No.: 24330	Expiration Date: 9/10/2026	Equipment Model: Guth 2100	
prescribed by the State Director of Tests pursu to 40 OAC 25-1-2.1, Approved breath-alcohol r	ared Spectroscopy utilizing the Intoxilyzer 8000 and w ant to 40 OAC 30-1-3. The above NIST traceable ref measurement equipment and reference solutions. Equ or Breath Alcohol Testers (72 FR 34747), published b	erence solution meets the requirements for use pursuant lipment used appears on the current or supplemented	
Signature Manual And	Date:	Date: 4/7/2025	
Name Christopher Pape	Permit No.: 2302	Permit No.: 23020	

RCVD 12/3/24

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24330 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on September 11, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0480% (w/vol) ethyl alcohol. The expiration date for this lot number is September 10, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06032102 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
01/15/2025 09:30

Test · ·	g/210L	Time
Air Blank	0.000	09:31
Cal Check	0.100	09:31
Air Blank	0.000	09:32
Cal Check	0.100	09:33
Air Blank	0.000	09:33
Cal Check	0.099	09:34
Air Blank	0.000	09:35
Cal Check	0.099	09:35
Air Blank	0.000	09:36
Cal Check	0.099	09:37
Air Blank	0.000	09:37
Cal Check	0.099	. 09:38
Air Blank	0.000	09:38
Cal Check	0.100	09:39
Air Blank	0.000	09:40
Cal Check	0.099	09:40
Air Blank	0.000	09:41
Cal Check	0.099	09:42
Air Blank	0.000	09:42
Cal Check	0.100	09:43
Air Blank	0.000	09:43
Cal Check Stats		
Average	0.0994	
Std Dev	0.0005	
Rel Std Dev(%)	0.5195	

Manufacturer: Guth	Manufactured Date: .	10/29/2024	Concentration: 0.100
Lot No.: 24390	Expiration Date: 10	/29/2026	Equipment Model: Guth 2100
prescribed by the State Director of Test to 40 OAC 25-1-2.1. Approved breath-a Conforming Products List of Calibrating	s pursuant to 40 OAC 30-1-3. The alcohol measurement equipment and	above NIST traceable refe d reference solutions. Equ 2 FR 34747), published by	as performed in accordance with the procedures erence solution meets the requirements for use pursuant ipment used appears on the current or supplemented y the National Highway Traffic Safety Administration.
Signature / August /		Date: 1/15/2025	
Name Christopher Pape		Permit No.: 23020	

RCVP 12/3/24

Certified Alcohol Reference Solution for Simulator

cp

Random Samples of Lot Number 24390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 30, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1205% (w/vol) ethyl alcohol. The expiration date for this lot number is October 29, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000

SN 80-002591

01/15/2025

09:45

Test	g/210L	Time
Air Blank	0.000	09:45
Cal Check	0.203	09:46
Air Blank	0.000	09:47
Cal Check	0.202	09:47
Air Blank	0.000	09:48
Cal Check	0.202	09:49
Air Blank	0.000	09:49
Cal Check	0.201	09:50
Air Blank	0.000	09:50
Cal Check	0.202	09:51
Air Blank	0.000	09:52
Cal Check	0.201	09:52
Air Blank	0.000	09:53
Cal Check	0.202	09:54
Air Blank	0.000	09:54
Cal Check	0.201	09:55
Air Blank	0.000	09:56
Cal Check	0.202	09:56
Air Blank	0.000	09:57
Cal Check	0.201	09:58
Air Blank	0.000	09:58
Cal Check Stats		
Average	0.2017	
Std Dev	0.0006	
Rel Std Dev(%)	0.3346	

lanufacturer: Guth	Manufactured Date: 2/13/2024	Concentration: 0.200
ot No.: 24080	Expiration Date: 2/13/2026	Equipment Model: Guth 2100
rescribed by the State Director of Tests 0 40 OAC 25-1-2.1. Approved breath-alc	a Infrared Spectroscopy utilizing the Intoxilyzer 8000 and versuant to 40 OAC 30-1-3. The above NIST traceable recohol measurement equipment and reference solutions. Equitis for Breath Alcohol Testers (72 FR 34747), published to Date:	ference solution meets the requirements for use pursuan julpment used appears on the current or supplemented by the National Highway Traffic Safety Administration.

12/3/24 07 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24150 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on April 17, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0241% (w/vol) ethyl alcohol. The expiration date for this lot number is April 16, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03122113 whose values are traceable to NIST.

Location Returned From	Unit Serial #	Date: 17/17
arvin Co 80-004885 4/11/25		
Unit checked in to BOT New Unit Checked out Checked out Remove old cylinder from	Old Cylinder #: 06023080A Unit Old cylinder disposed in inv n unit (After Lab Evaluation) Old c	ventory
Check in Procedure Connect unit to ethernet connection 172.019.004. 137 Download any stored records 2 Change unit location to Serial number Update Unit Date and time Previous time: 4/14/2035 1319 New Time: 7 Me OK Execute a new Forms load on the unit	Notes from the Field ✓ Routine Swap	ssues in the Field
Lab Evaluation Breath Hose Hold Vacuum?: Test Ch	namber Hold Vacuum?:	
DVM 3μm 13118 9μm 14095 Peripheral Tests Chamber Temp 47.0 Breath Hose Temp 42.0	Snap Test ok?	
Speaker (F5) Power LED Green (F6) Gas Regulator (F11) Pump (F10)	Red (F7) Fan (F9)	
Atmospheric Sensor Pressure 972 Correction Factor 1.0	Flow Sensor 200	Is regulator close to Tank Sensor Value?
Sent to CMI (If needed) Reason Sent Date Sent Invoice Sent Date Returned		
Repairs needed:	Repairs Made:	
	Ready for Bench Check? ✓	
Inventory Process □ Bench Check Techn □ Linked Cylinder to Unit □ Mock Subject Test Oralism □ Changed Unit Status □ Certificate of Calibration □ Bench Check Report Completed □ Service Overview Completed	Completed CMI Workord Cation and Operation CMI Certificat	te of Calibration