



**STATE OF OKLAHOMA  
BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE**

Dr. Kenneth E. Blick, Ph.D.  
Chairman

Administrative Office at:  
Board of Tests for Alcohol and Drug Influence  
Post Office Box 36307  
Oklahoma City, Oklahoma 73136-2307

J. Robert Blakeburn  
State Director of Tests  
for Alcohol and Drug Influence

---

Emmett E. Goodson  
Administrative Assistant


**April 15, 2008**

**ACTION NO. 08-03 BY THE STATE DIRECTOR OF TESTS FOR ALCOHOL AND DRUG  
INFLUENCE**

Under authority of 47 O.S. Supp.2000, Section 751-761 and OAC 40:1-1-3. and 40:30-1-3. Breath-alcohol analysis with the Intoxilyzer Model 8000; promulgated by the Board of Tests for Alcohol and Drug Influence, the following official **ACTION** is hereby taken by the State Director of Tests for Alcohol and Drug Influence.

**LOCATION REQUIREMENTS FOR OPERATING THE MODEL 8000 INTOXILYZER**

- (1) The State Director of Tests for Alcohol and Drug Influence hereby prescribes and approves the following requirements for locations utilizing the CMI Intoxilyzer 8000 Breath instrument.
- (2) The CMI Intoxilyzer 8000 and printer should be placed in a location which is:
  - (a) **Clean.** The area should not be used for storage of, or commonly subject to vapors of, volatile materials.
  - (b) **Secure.** While not in use the instrument and printer should be secure from unauthorized personnel.
- (3) The preceding requirements are prospective in nature and are intended to prevent damage to instruments and delays in testing. If an instrument prints an affidavit after a test has been administered, it indicates the requirements were met at the time the test was administered.
- (4) This **ACTION** by the State Director of Tests for Alcohol and Drug Influence shall take effect on and after April 15, 2008.

  
J. Robert Blakeburn  
State Director of Tests  
for Alcohol and Drug Influence